

## Determination of estriol in serum during pregnancy

A-22.2

### Key words

Instrumental HPTLC - quantitative analysis - prechromatographic derivatization - densitometry (fluorescence) - clinico-chemical diagnostics - steroids - estriol

### Scope

Estriol is present in serum mainly in form of its conjugates (16-glucuronide, 3-glucuronide, 3-sulfate and 3-sulfate-16-glucuronide).

Estriol is liberated from its conjugates by  $\beta$ -glucuronidase, extracted by solid phase, and derivatized with dansyl chloride. Following chromatography on silica gel with toluene - dioxane - methanol 8:2:1, it is determined quantitatively by fluorescence measurement.

The limit of quantification of this procedure is 2  $\mu\text{g/L}$  estriol in serum. The results correlate outstandingly well with those obtained by radioimmunoassay (RIA).

### Literature

W. Funk, Fresenius Z. Anal. Chem. **318**, 206-219 (1984).

F. Arndt, Thesis, Giessen University, Faculty of Human Medicine, 1983.

#### Advantages of using HPTLC for this analytical task

- Substantially greater cost effectiveness than RIA, particularly if the daily (weekly) number of such analyses is only small.

## Chemicals

|  |                                   |
|--|-----------------------------------|
| β-Glucuronidase (Merck) 40 U/mL          | Proline (3.5 mg/mL water bidist.) |
| Sodium hydroxide 1N                      | Methanol                          |
| ClinElut (Analytichem Int. USA No. 1001) | n-Hexane                          |
| Dichloromethane                          | Toluene                           |
| Isopropanol                              | Dioxane                           |
| Acetone                                  | Paraffin oil                      |
| Dansyl chloride                          |                                   |

Standard: Estriol

## Sample preparation and derivatization

a) Enzymatic cleavage and extraction:

- Add 20 μL β-glucuronidase to 500 μL serum, homogenize and incubate the mixture for 2 h at 37°C.
- Adjust to pH 8.5 with 1N sodium hydroxide, pour onto ClinElut cartridges, and allow to adsorb for 3 min.
- Elute with 4 mL dichloromethane - isopropanol 9:1, repeat it after 4 min.
- Evaporate the eluate to dryness in a stream of nitrogen.
- Dissolve residue in 1.5 mL acetone, transfer to an Eppendorf reaction vessel, and evaporate again to dryness.
- Dissolve residue in 100 μL acetone.

b) Derivatization:

- Mix the residue dissolved in 100 μL acetone with 20 μL 0.2% dansyl chloride solution and 10 μL 1N sodium hydroxide, shake the mixture for 15 s.
- Add 50 μL acetone, shake, and incubate for 60 min in the dark at room temperature.
- Add 20 μL proline solution, swirl and incubate again for 15 min in the dark at room temperature. (The excess dansyl chloride dissolves).
- Use supernatant phase for application onto the HPTLC plate.  
(10 μL = 25 μL human serum).

## Standard solution

Dissolve 10 mg estriol in 100 mL acetone, dilute 1:100 with acetone and derivatize in the same way as the samples.

## Layer

HPTLC plates silica gel Merck 60 F<sub>254</sub>, 20x10 cm, prewashed with methanol.

### Sample application

With CAMAG Linomat as 7 mm bands, distance between tracks 3 mm, distance from left edge 15 mm, distance from lower edge 7 mm, delivery rate 5 s/ $\mu$ L = 18 applications per plate.

Application pattern:

|    |    |    |    |    |    |    |    |    |    |                  |
|----|----|----|----|----|----|----|----|----|----|------------------|
| S1 | U1 | S2 | U2 | S3 | U3 | S4 | U4 | S1 | U1 | ...              |
| 1  | 10 | 2  | 10 | 4  | 10 | 8  | 10 | 1  | 10 | ... $\mu$ L/Band |

S1- S4 = standard in different concentrations; U1-U4 = unknowns

Alternatively the Automatic TLC Sampler can be used.

### Chromatography

In CAMAG Horizontal Developing Chamber 20x10 cm with toluene - dioxane - methanol 8:2:1 with chamber saturation, migration distance 5 cm, running time 13 min,  $R_f$  about 0.36.

While still wet, dip the plate with the CAMAG Chromatogram Immersion Device for 1 s in a solution of paraffin oil - n-hexane 5:1. The fluorescence intensity is increased tenfold and stabilized for several hours.

### Densitometric evaluation

With CAMAG TLC Scanner and CATS evaluation software; scanning by fluorescence at 313/>400 nm.

The limit of quantification for estriol is 100 pg/zone, corresponding to a serum concentration of about 2  $\mu$ g/L.

### Discussion

Table 1: Estriol level in the serum ( $\mu$ g/L)

| Total estriol<br>(conjugated and unconjugated) | unconjugated estriol |        |
|--|----------------------|--------|
| normal cycle                                   | 0.25 - 0.37          |        |
| week 25 of pregnancy                           | 50                   | 2 - 6  |
| week 40 of pregnancy                           | 200                  | 7 - 25 |




Fig. 1 Calibration curve for estriol




Fig. 2 Comparison with radioimmunoassay

Results obtained by quantitative HPTLC determination of 10 different sera compared with those found by RIA. Orthogonal regression yields a correlation coefficient of  $r = 0.9778$  and confirms the good correlation of the two analytical procedures.