

## Determination of lactose, saccharose and fructose/glucose

A-25.2

### Key words

Instrumental HPTLC - quantitative analysis - densitometry (absorbance) - post chromatographic derivatization - process control - food analysis - carbohydrates - biosynthesis-monitoring - antibiotics

### Scope

The method was developed to monitor the biosynthesis of antibiotics. Fermentation broth (culture medium) is diluted, filtered and chromatographed without further processing. The plate is then derivatized with aniline-diphenylamine-phosphoric acid reagent, and evaluated densitometrically.

The accuracy attainable is 2-3 % (CV at n = 10 on different plates).

### Literature

F. Kreuzig, J. Liquid Chromatog. **6**, 1227-1236 (1983)

#### **Advantages of using planar chromatography for this analytical task**

- Extremely easy sample preparation
- Low running costs
- Speed and reliability of method
- Elimination of interferences that can be caused by enzyme inhibitors in enzymatic methods

## Chemicals

Methanol	Standards:
Water	Lactose
1-Butanol	Saccharose
Acetic acid	Fructose
Aniline	Glucose
Phosphoric acid 85%	

## Sample preparation

The culture medium (fermentation broth) is diluted with methanol - water 1:1 until the concentration (total sugar) is about 0.5%. The solution is filtered and chromatographed immediately.

## Standard solutions

Dissolve 10 mg each of lactose, saccharose, fructose or glucose\* to 10 mL in methanol - water 1:1 and dilute this solution 1:2 with methanol - water 1:1 (500 ng/μL).

Application scheme:

U1	S1	U2	U3	S2	U4	U5	S3	U6	U1	S1	U2	...
	200			600			1000			200		μL/band
	100			300			500			100		ng carbohydrate each

## Layer

HPTLC plates Merck silica gel 60 F<sub>254</sub>, 20 x 10 cm

## Sample application

Bandwise with CAMAG Automatic TLC Sampler III, band length 4 mm, distance between tracks 6 mm, distance from side 15 mm, distance from lower edge 10 mm = 17 applications.

## Chromatography

In CAMAG Twin-Trough chamber with 1-butanol - acetic acid -water 80:100:15 with chamber saturation, migration distance 50 mm, running time about 30 min.

After chromatography, the plate is dried for 20 min at 115°C. Then the chromatogram can be stored. The following derivatization is carried out when densitometric evaluation is feasible within one hour.

---

\* In this method fructose and glucose are not separated. If they are to be determined separately, application method A-07.3 (NH<sub>2</sub>-bonded silica) must be used.

## Derivatization

By dipping for 3 s with CAMAG Chromatogram Immersion Device in diphenylamine reagent (3 mL aniline + 3 g diphenyl amine + 15 mL conc. phosphoric acid made up to 150 mL with methanol) followed by heating at about 115°C for about 15 min. The colors are as follows:

Lactose: greyish blue


Saccharose: brownish green

Fructose: rust brown

(Glucose lies between saccharose and fructose.)

## Densitometric evaluation

With CAMAG TLC Scanner and CATS evaluation software; scanning by absorbance at 620 nm.



Densitogram of 1 lactose, 2 saccharose, 3 fructose