




Cartridger® C-670

High efficiency cartridge packing

The patented Cartridger® C-670 allows quickly and simply packing high performance cartridges for reliable flash chromatography separations. Together with a substantial cost reduction, it provides the possibility to fill the cartridge with the silica quality of your choice.

Unique

Proprietary method for high efficiency packing



Cost saving

Economical solution for your purifications



Fast and safe

Fastest process and highest user protection

C-670: Your most important benefits



Unique

- Achieve optimal separation with high quality cartridges
- Benefit from the unique packing process based on the fluidization of the silica
- Prepare your specific cartridge: choose your own silica and select the cartridge format (4 to 90 g)



Cost saving

- Save cost compared to prepacked cartridges
- No more column storage
- Reduce solvent consumption compared to open glass columns



Fast and safe

- Pack a column in less than one minute
- Use any silica based material with particle size from 30 to 200 μm
- Minimized user exposure

Complete your portfolio



Pump Module
C-601 / C-605
Pump



Fraction collector
C-660
Fraction collection



UV-Vis Detector
C-640
UV-Vis detection



Glass columns
Columns for scale up

