



Pump Module C-601 / C-605

Reliable solvent delivery

The pump modules C-601 and C-605 are the core components of the BUCHI Sepacore® Chromatography systems. Whatever the application in isocratic or gradient mode is, the modules provide a reproducible and pulse free flow rate.



Universal

Pump any solvent at low or medium pressure



Small and powerful

Smallest footprint with broad flow rate range and gradient capability



Flexible

Upgradeable according to your application and budget



C-601 / C-605: Your most important benefits



Universal

- Free choice of solvents thanks to inert pump head material
- Reliable flow rate up to 50 bar
- Multiple control modes according to the application



Small and powerful

- Smallest footprint with flow rate up to 250 mL/min
- Pulse free solvent delivery using a 3-piston pump head
- Plug and play portable unit



Flexible

- Configuration adjusted to your needs and budget
- Easily upgradeable: isocratic to gradient or manual to computer control

Complete your portfolio



Fraction collector
C-660
Fraction collection



Sepacore® Flash
cartridges
Prepacked columns



UV-Vis Detector
C-640
UV-Vis detection



SepacoreControl
Software
Full control software

