



Distillation Unit K-350 / K-355

Entry level Kjeldahl analysis and beyond

Perform Kjeldahl nitrogen or protein determinations compliant to officially approved methods in the most cost effective way. Expand your range of applications with the additional features of the K-355.



Cost-effective

Rapid steam distillation
within your budget



Easy to use

Simple programming and
straightforward operation



Wide-range use

Comply with approved Kjeldahl
methods and extend your
application range

K-350 / K-355: Your most important benefits



Cost-effective

- Entry level instrument for Kjeldahl and other steam distillations
- Eco-friendly cooling water control to reduce water consumption
- Plastic splash protector for extended lifetime



Easy to use

- Simple programming of distillation parameters
- Straightforward operation and calibration of reagent pumps
- Multiple method memory for routine analysis with the K-355



Wide-range use

- Acid resistant pump for additional steam distillation applications
- Reduction of steam power for gentle distillation
- Choice of glass or plastic splash protector

Complete your portfolio



SpeedDigester
K-439 / K-425 / K-436
IR digestion



Extraction Unit
E-816 SOX & HE
Dedicated fat extraction



NIRMaster™
IP54 / Pro IP65
FT-NIR spectroscopy



Recirculating Chiller
F-108
The efficient way of
cooling

