



THERMOSTAT

-10+50°C / 0,1°C
SILENT



- 1. Proportional derivative electronic temperature control from -10 to +50 °C with accuracy ± 0.1 °C**
- 2. Easy temperature setting**
- 3. Useful space of STAINLESS STEEL TANK 30x35x19 cm (9 bottles)**
- 4. Cooling coil stainless steel**

5. Power supply 230V 1600W.


6. Dimensions: (W X D X H) cm 32 x 58 x 83

7. Weight: 39 Kg

Pivoting wheels as standard.

SPECIAL EXECUTIONS on REQUEST


Accessory on request: pump for external circulation of the liquid to drain the tank.





LARGE USEFUL SPACE
IN THE TANK ABOUT 9
BOTTLES

The characteristics of this documentation are not binding and may change without notice.
All our instruments are in compliance with Directive 2004/108 on electromagnetic compatibility.


customer care
info@gibertini.com
Tel. (+39) 023541434

ISO 9001:2008

ACCREDIA
L'INTERNAZIONALE ITALIANA
LAT N° 094

GIBERTINI
ELETTRONICA s.r.l.
Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (+39) 02 3541434
Fax (+39) 02 3541438
www.gibertini.com
sales@gibertini.com