

# SPECIAL OFFER



ADDITIONAL 10% DISCOUNT  
FOR UNIVERSITIES



## TINYCLAVE

steel vessel 10 ml  
up to 100 bar  
-20 to +300 °C

**CHF 2'310**

## MINICLAVE

steel vessel 100 ml  
up to 100 bar  
-20 to +300 °C

**CHF 2'780**

- Interchangeable reaction vessels available in glass, stainless steel (insert PTFE) and Hastelloy which are all compatible with the same reactor lid. Used for chemical reactions up to 10 bar in glass and up to 100 bar in metal reaction vessels.
- Volume range tinyclaves: 10 to 25 ml  
Volume range miniclaves: 100 to 300 ml